



# Generating a Single Cell Suspension from Adherent Cells

## OVERVIEW

Specific cell culture procedures, and products to allow for Cancer Stem Cell isolation and propagation as tumor spheres.

### DISRUPTION OF ADHERENT CELL PROCEDURE

Starting from a confluent (90%) culture in a T75 tissue culture flask

1. Remove existing medium, and rinse adherent cells with 5 mL of sterile Hanks Balanced Salt Solution (HBSS), or Phosphate Buffered Saline (PBS); discard wash (this will remove trace amounts of fetal bovine serum (FBS), which will normally may inhibit the action of trypsin).
2. Add 5 mL of pre-warmed (37°C) Trypsin:EDTA, and place the flask in a 37C incubator for 2 min.
3. Following incubation cells should be loosely attached to the flask surface, and cells are disrupted by repeated pipetting (using a 5 mL sterile pipette).

*The Trypsin: EDTA incubation time is critical to avoid excess clumping and cell death. While 2 min will work for many cell types/lines, this may need to be optimized based on the specific cell type/line being used.*

4. Following disruption, immediately add 15 mL of DMEM media, containing 10% FBS, to the trypsinized cells, transfer to a sterile 50 mL polypropylene conical centrifuge tube, and pellet cells by centrifugation at 1000 RPM for 5 min.
5. Resuspend cells in 10 mL of Cancer Stem Premium™ (cat. # 20101, ProMab Biotechnologies) media.
6. The cell suspension is passed through a sterile 20 µM cell strainer to remove any cell clumps remaining in the cell suspension.
7. Remove a 10 µL aliquot for cell counting.

#### Additional Products and Services:

-  [Mouse Monoclonal Antibody](#)
-  [Rat Monoclonal Antibody](#)
-  [Rabbit Monoclonal Antibody](#)
-  [Human Monoclonal Antibody](#)
-  [Polyclonal Antibody](#)
-  [Antibody Sequencing](#)
-  [Hybridoma Sequencing](#)
-  [CAR T-cells](#)
-  [Lentivirus production](#)
-  [Cancer Stem Cells](#)
-  [Specialty Cell Culture Media](#)
-  [T-cell Expansion beads](#)

Ask about our full line of CRO services to provide supplemental assistance or the entire support necessary to complete your project on time and with the data you need to move forward.