



Dengue NS1, ST4 Recombinant Protein

Product Information

Catalog#	Pr20483
Name	Dengue NS1, ST4 Recombinant Protein

Order Information

Size	1000µg	500µg
Price	\$1200	\$600

Specification

Aliases	
Storage	Five years frozen. One month in solution at room temperature.
Formulation	(1mg/1.3ml) 100mM Tris HCl pH-8, 8M Urea & 200mM NaCl.
Category	Viral Antigens
Purity	Protein is >90% pure as determined by 10% PAGE (coomassie staining).
Source	













Description

The E.Coli derived recombinant protein contains the NS1 Dengue Virus full length Type-4 immunodominant regions. The dengue protein is fused to 6xHis tag.

Sequence

Physical

Additional Products and Services:

-  [Mouse Monoclonal Antibody](#)
-  [Rat Monoclonal Antibody](#)
-  [Rabbit Monoclonal Antibody](#)
-  [Human Monoclonal Antibody](#)
-  [Polyclonal Antibody](#)
-  [Antibody Sequencing](#)
-  [Hybridoma Sequencing](#)
-  [CAR T-cells](#)
-  [Lentivirus production](#)
-  [Cancer Stem Cells](#)
-  [Specialty Cell Culture Media](#)
-  [T-cell Expansion beads](#)

Ask about our full line of CRO services to provide supplemental assistance or the entire support necessary to complete your project on time and with the data you need to move forward.