



YWHAQ Human Recombinant Protein

Product Information

Catalog#	Pr22507
Name	YWHAQ Human Recombinant Protein

Order Information

Size	1mg	20µg
Price	\$3600	\$130

Specification

Aliases	14-3-3 theta, 14-3-3 tau, 14-3-3 T-cell, HS1, YWHAQ, 1C5, 14-3-3 Tau, Tyr-3/Trp- 5 Monoxygenase Activation Protein Theta.
Storage	YWHAQ Human Recombinant although stable at 4°C for 1 week, should be stored desiccated below -18°C. Please prevent freeze thaw cycles.
Formulation	YWHAQ solution containing 20mM Tris 7.5.
Category	Enzymes
Purity	Greater than 95.0% as determined by(a) Analysis by RP-HPLC.(b) Analysis by SDS-PAGE.
Source	Escherichia Coli.

Description

YWHAQ Human Recombinant produced in E.Coli is a single, non-glycosylated, polypeptide chain containing 245 amino acids (1-245) and having a molecular mass of 27 kDa. YWHAQ is purified by proprietary chromatographic techniques

Sequence

MEKTELIQKA KLAEQAERYD DMATCMKAVT EQGAELSNEE
RNLLSVAYKN VVGRRSAWRVISSIEQKTD TSDKKLQLIK
DYREKVESEL RSICTTVLEL LDKYLIANAT
NPESKVFYLMKGDYFRYLA EVACGDDRKQ TIDNSQGAYQ
EAFDISKKEM QPTHPIRLGL ALNFSVFYIEILNNPELACT
LAKTAFDEAI AELDTLNE

Physical

Sterile filtered colorless solution.

Additional Products and Services:

-  [Mouse Monoclonal Antibody](#)
-  [Rat Monoclonal Antibody](#)
-  [Rabbit Monoclonal Antibody](#)
-  [Human Monoclonal Antibody](#)
-  [Polyclonal Antibody](#)
-  [Antibody Sequencing](#)
-  [Hybridoma Sequencing](#)
-  [CAR T-cells](#)
-  [Lentivirus production](#)
-  [Cancer Stem Cells](#)
-  [Specialty Cell Culture Media](#)
-  [T-cell Expansion beads](#)

Ask about our full line of CRO services to provide supplemental assistance or the entire support necessary to complete your project on time and with the data you need to move forward.