

## PM-CAR2013 NON-TRANSUCED MACROPHAGES

Ready-to-use Macrophage cells

Primary macrophages can be isolated from PBMC, monocytes using Percoll-Ficoll purification, monocyte specific markers and differentiated with PMA agent. The macrophages can be used for different functional tests: phagocytosis, differentiation, polarization, etc

### Data




Figure 1. Primary macrophages induce Jurkat cell phagocytosis in the presence of CD47 antibody.

#### Products and Services

- Mouse Monoclonal Antibody
- Rat Monoclonal Antibody
- Human Antibody
- Hybridoma Sequencing
- Polyclonal Antibody